

**Department of Physics and Astronomy**  
**University of Heidelberg**

Bachelor Thesis in Physics  
submitted by

**Robin Eberhard**  
born in Aalen, Germany

handed in on  
**March 28, 2016**



# **Characterization of a multispecies imaging system**

This Bachelor Thesis has been carried out by Robin Eberhard at the  
**Physikalisches Institut in Heidelberg**  
under the supervision of  
**Prof. Dr. Matthias Weidemüller**



# Characterization of a multispecies imaging system

Robin Eberhard

**Abstract** Lorem ipsum dolor sit amet, consectetur adipiscing elit. Ut purus elit, vestibulum ut, placerat ac, adipiscing vitae, felis. Curabitur dictum gravida mauris. Nam arcu libero, nonummy eget, consectetur id, vulputate a, magna. Donec vehicula augue eu neque. Pellentesque habitant morbi tristique senectus et netus et malesuada fames ac turpis egestas. Mauris ut leo. Cras viverra metus rhoncus sem. Nulla et lectus vestibulum urna fringilla ultrices. Phasellus eu tellus sit amet tortor gravida placerat. Integer sapien est, iaculis in, pretium quis, viverra ac, nunc. Praesent eget sem vel leo ultrices bibendum. Aenean faucibus. Morbi dolor nulla, malesuada eu, pulvinar at, mollis ac, nulla. Curabitur auctor semper nulla. Donec varius orci eget risus. Duis nibh mi, congue eu, accumsan eleifend, sagittis quis, diam. Duis eget orci sit amet orci dignissim rutrum.

**Zusammenfassung** Lorem ipsum dolor sit amet, consectetur adipiscing elit. Ut purus elit, vestibulum ut, placerat ac, adipiscing vitae, felis. Curabitur dictum gravida mauris. Nam arcu libero, nonummy eget, consectetur id, vulputate a, magna. Donec vehicula augue eu neque. Pellentesque habitant morbi tristique senectus et netus et malesuada fames ac turpis egestas. Mauris ut leo. Cras viverra metus rhoncus sem. Nulla et lectus vestibulum urna fringilla ultrices. Phasellus eu tellus sit amet tortor gravida placerat. Integer sapien est, iaculis in, pretium quis, viverra ac, nunc. Praesent eget sem vel leo ultrices bibendum. Aenean faucibus. Morbi dolor nulla, malesuada eu, pulvinar at, mollis ac, nulla. Curabitur auctor semper nulla. Donec varius orci eget risus. Duis nibh mi, congue eu, accumsan eleifend, sagittis quis, diam. Duis eget orci sit amet orci dignissim rutrum.



# Contents

<b>1. Introduction</b>	<b>1</b>
<b>2. Setup for high resolution imaging</b>	<b>3</b>
2.1. Experimental requirements . . . . .	3
2.2. Camera for double species imaging . . . . .	4
2.2.1. Basics of CCD cameras . . . . .	4
2.2.2. Comparison with the present setup . . . . .	6
2.2.3. Dark current . . . . .	7
2.2.4. Readout noise . . . . .	8
2.2.5. Quantum efficiency . . . . .	10
2.3. Mechanical shutter . . . . .	10
2.3.1. Electronic and mechanical setup . . . . .	11
2.3.2. Dynamical properties . . . . .	12
2.4. Mask for the CCD sensor . . . . .	15
2.4.1. Fast kinetics mode . . . . .	15
2.4.2. Frequency response of a slit . . . . .	17
2.4.3. Optimization of the masking setup . . . . .	20
<b>3. Thermometry of an ultracold ideal fermi gas</b>	<b>23</b>
3.1. Absorption imaging . . . . .	23
3.2. Density distributions of ideal Fermi gases . . . . .	25
3.3. Finding properties of the Fermi gas . . . . .	26
<b>4. Conclusion and outlook</b>	<b>29</b>
<b>A. Acquisition sequence</b>	<b>31</b>
<b>B. Shutter circuit</b>	<b>33</b>
<b>C. Setup of the custom slit</b>	<b>35</b>





# 1. Introduction

Conducting experiments, measuring attributes and taking data are quantities that go hand in hand in scientific fields. Already in the early 1900s, imaging systems were chosen as data acquisition in astronomy. A complete optical system allowed for magnification of planets, our sun or nebulae. This trend started with the use of analogous photo plates [1], while nowadays so-called charge-coupled devices (CCD) or complementary metal-oxide-semiconductors (CMOS) are used as they can take the data digitally.

This also found application in the field of ultracold quantum physics — an experimental field that is around since the 1970s. Hereby, atoms are prepared in traps in order to reduce their velocities. The atoms can then be detected by illuminating them with an imaging beam, a laser beam that is resonant with the atoms. The atoms will cast a shadow, which can be focused onto the chip, where the result is digitalized and can be analyzed. This process offers some challenges. At low temperatures, the trap width can not be arbitrarily big. This means, that only fewer atoms can be imaged, therefore a high resolution imaging is required in order to distinguish them from the background noise. But the resolution of the system is not only limited by the chip and pixel size of the camera, it also needs to be refined on optical elements such as lenses and mirrors.

The experiment at this group uses a mixture of  $^6\text{Li}$  and  $^{133}\text{Cs}$ , in order to measure for example the Efimov effect, ideal Fermi gases or Polarons. An important part in this setup is the acquisition, where the species need to be acquired separately. A new system now allows to image both  $^6\text{Li}$  and  $^{133}\text{Cs}$  only after a short delay, instead of taking the images from different axes or clouds.

This thesis builds upon the optical system, that made the high resolution possible and will introduce a new camera into the experiment, which has new features, that allow for noise reduction and an improvement in acquisition timings.



## 2. Setup for high resolution imaging

In order to understand the nature of atoms in more detail, experiments are conducted where they are trapped and cooled. Under these circumstances, interactions between atoms can be tuned using Feshbach resonances. This allows to build up for example Bose-Einstein condensates (BEC) or Cooper pairs. The complete setup has been explained in more detail before [2]. The following sections will focus on the imaging in order to extract these attributes and the measures it takes to reduce noise sources on the camera.

### 2.1. Experimental requirements

In this double species experiment, fermionic  $^6\text{Li}$  and bosonic  $^{133}\text{Cs}$  are trapped and cooled inside a vacuum chamber. The atoms are emitted from an oven into a Zeeman slower, which transfers them into a magneto-optical trap. The MOT will further cool the atoms up to the doppler limit ( $140\text{ }\mu\text{K}$ ). To reach lower temperatures, the atoms will be trapped inside an optical dipole trap. Forced evaporative cooling will release the fastest atoms, until temperatures of a few nano Kelvin are reached.




Figure 2.1.: **Imaging path.** The camera is mounted above the vacuum chamber as visible in **a.** (gravity would be in negative  $z$  direction). The image of the atoms is collimated with the first lens and passing the optics in **b.**, until it is refocused in the chip on the camera.

In order to image the atoms, an imaging beam is pointing along the  $z$ -axis according to Figure 2.1. An achromatic doublet lens focuses the imaging beam onto the CCD camera, which is located on top of the vacuum chamber.

In order to measure even at low atom numbers ( $n = 1000$ ), the setup was refined allowing for high resolutions. The first lens, which collimates the image of the atoms has a low focal length ( $f_1$ ). The image is refocused into the camera by the second lens. Since the imaging uses two separate frequencies to image both Lithium and Caesium, the chromatic shift introduced by the lenses will result in two different focal points for the images. This can be compensated by a high focal length, which therefore allows to have the camera in the Rayleigh range of both laser beams.

This setup then allows to image for example Bose-Einstein gases, polarized Fermi gases or superfluids with high accuracy.

## 2.2. Camera for double species imaging

The camera used in order to find atoms, the Andor iKon M, is a charge-coupled device (CCD), which offers an improvement to the old setup. Using this camera, noise sources such as Readout noise and Dark noise can be reduced, which are explained in the following sections.

### 2.2.1. Basics of CCD cameras

A camera operates by means of converting photons first into electrons then into voltage, which is finally read out as data [3] [4]. Each conversion process can add noise to the final image, which needs to be minimized in order to acquire accurate data. The most common noise sources are covered in Chapter 2.2.

The photons are collected on an array of semiconductor photo diodes, called the pixels, where ideally the spacing between the pixels is zero to gain maximum accuracy. The resolution is then dependent on the pixel size, which is for scientific cameras usually between  $10\text{ }\mu\text{m}$  and  $20\text{ }\mu\text{m}$  per pixel. Bigger pixels means higher photon sensitivity but usually lower resolution.

To create a digital image, the charge from the pixels have to be shifted one-by-one into the analog digital converter (ADC). This is done, by vertically shifting them into the readout

get back to  
stephan. is  
this correct?

focal length




Figure 2.2.: **Schematic design of a CCD array and its readout.** The pixels are arranged in the pixel array. During readout they are shifted into the readout register and then to the side into the analog to digital converter

register and then horizontal into the ADC, where the charges are multiplied and converted to digital data.

The shifting is done by storing the charges after collecting them, where each storage can be seen as an electronic potential well. In order to shift the charges and prevent overlapping, three spatially separated regions with potentials  $U_1$ ,  $U_2$  and  $U_3$  are required. Figure 2.3 indicates the systematics behind the shifting.




Figure 2.3.: **Shifting charges in a CCD detector.** To shift charges from Pixel A to Pixel B, the three potentials at the regions  $U_1$ ,  $U_2$  and  $U_3$  in each pixel have to be set accordingly to allow the charge flow without overlapping each other. Each row represents a single step.

Moving electrons to the next region is a three-step process. The charges are first only

present in the region  $U_2$ , while the potentials  $U_1$  and  $U_3$  are kept high. They are then distributed across  $U_2$  and  $U_3$  by setting them both low. At last  $U_2$  will be set high such that the charge is now fully in  $U_3$ . This is repeated until the charge has been finally shifted from pixel A to B, which in total takes seven steps.

### 2.2.2. Comparison with the present setup

For the imaging of small atomic clouds, it is very important to have cameras with optimal noise and maximal readout speed. The new setup improves on both attributes. The dark noise can be significantly reduced by cooling the chip down to more than 70 °C. A new readout mode, called fast kinetics, makes it possible to acquire all images before reading out. This significantly improves the speed at which images can be taken.

The readout speed is highly important in our setup. Since absorption imaging is the technique of choice to measure atom attributes, three images need to be taken each sequence, being the absorption, division and background image. The old setup used a Guppy-38B [2] camera, which has a frame rate of 30fps. This meant the acquisition was finished after 100 ms. The Andor camera, on the other hand, can take images quickly without the need to read out in-between. At the fastest shift speed, the acquisition is finished after 1.632 ms, improving the speed by a factor of more than 60.

The quantum efficiency is explained in Chapter 2.2.5 and is also highly important, as it describes how many photons are detected on a camera and is directly connected to the sensitivity. A higher quantum efficiency therefore means better results. They compare as  $QE_{Li} = 79\%$ ,  $QE_{Cs} = 77.3\%$  for the Andor and  $QE_{Li} = 35\%$ ,  $QE_{Cs} = 10\%$  for the Guppy camera, therefore the sensitivity of the new camera is significantly higher.

When comparing the resolution, the chip size also has to be considered. Since higher resolutions seem to be preferable at first, it also means that for the same pixel sizes, the photon sensitivities will decrease. The pixel sizes compare to  $8.4\mu\text{m} * 9.8\mu\text{m}$  to  $13\mu\text{m} * 13\mu\text{m}$  from the old versus the new setup respectively, while the resolutions are  $768 \times 492$  and  $1024 \times 1024$  making a larger magnification now possible.

The Guppy camera is a lot smaller than the Andor iKon M ( $48.2\text{ mm} \times 30\text{ mm} \times 30\text{ mm}$  vs  $204.2\text{ mm} \times 105\text{ mm} \times 107\text{ mm}$ ), therefore making an implementation on a full experimental table easier. To implement such a complicated camera system that the Andor is, a lot of preparation was made in the thesis of Carmen Renner [5].

Nevertheless, despite its size, the new camera offers the ability to image both Lithium and Caesium species at once, while two Guppy cameras were needed beforehand, which also meant placing them on different imaging axes.

### 2.2.3. Dark current

A common noise source that is apparent in all CCD cameras, is the so-called dark current [6]. It originates from the thermal excitation of electrons in individual pixels. Since they are made of semiconductors, once in a while, an electron can pass the potential between valence and conduction band simply due to their thermal energy. Thus, excess electrons accumulate, that contribute to the background signal and introduce additional noise [7].

The dark current has a strong temperature dependence

$$I_{dark}(T) \propto T^{\frac{3}{2}} \exp(-E_g/2k_B T), \quad (2.1)$$

while  $E_g$  is the band gap, that separates the valence from the conduction band in the semiconductor,  $T$  the temperature and  $k_B$  the Boltzmann constant. Therefore in order to reduce dark current noise, the temperature of the chip can be reduced, which decreases the thermal energy of the electrons.

This has been measured and verified in Figure 2.4, using the built-in peltier element in the Andor camera to control the temperature of the CCD chip. For a long exposure time, dark current accumulates on all pixels. The counts are measured for several temperature settings and the counts  $C$  are converted to electrons per pixel per second ( $I_i$ ) [6] for each pixel  $i$  as

$$I_i = \frac{C}{G t_{exp}} \quad (2.2)$$

with the Gain  $G$  of the camera, and the exposure time  $t_{exp}$ . The readout time of the pixels can be neglected, as they are significantly lower than the exposure time (see Chapter 2.4.1). The data should then follow the theory in Equation 2.1, although deviations are visible in the low temperature regime, which result from the fan not being able to divert the heat from the chip.

As it is important to minimize the noise source, the chip is cooled as low as possible. This will reduce the possibility of thermal electrons creating counts in the detector, so that even at fewer atom numbers, they are still distinguishable from the noise.




Figure 2.4.: **Dark noise accumulation on the chip.** For a long exposure time of 100 s, the dark current was measured for several temperatures. Gain in this measurement was minimal (0.215 [5]). Equation 2.1 was fitted, leaving the amplitude  $A$  and band gap  $E_g$  free, which resulted in  $A = 1.1 \cdot 10^8 \text{ W m}^{-1} \text{ K}^{-1}$  and  $E_g = 0.96 \text{ eV}$ .

## 2.2.4. Readout noise

As described in Chapter 2.2.1, pixels are shifted in order to be read out by the ADC. Moving charges from pixel to pixel causes noise that accumulates over each iteration. Together with the noise which the ADC introduces, this is then called the readout noise. The shifting noise together with the dark noise is visible as a gradient () since each shift adds new charges due to excitations in the semiconductors. In order to characterize this, one can take the variance of the image, which should be zero for no noise.

In Figure 2.5 it has been shown, that slow readout speeds reduces the noise, although the vertical shift speeds do not seem to influence the noise at all. The measurement was carried out by taking images at low exposure time and low temperature, so that the main noise source would be the readout and not the dark current. Therefore the variance in the plot shows the noise from readout which goes to zero as the horizontal speed is decreased.

The actual readout of the chip can be done at slow speeds, so that the noise from horizontally shifting is low. As the vertical shifts do not add significant noise, it can safely be set to fast speeds, which is also favorable, as  $^6\text{Li}$  and  $^{133}\text{Cs}$  are imaged back to back as fast as possible, so that the cloud dynamics are still approximately the same.

Shifting charges means increasing and decreasing potential wells, so that they can move from one into the other. If this is done fast enough (non-adiabatic), electrons in the pixel are more likely to get excited from the valence to the conduction band from the semiconductor,

figure with  
shift gradi-  
ent



therefore adding noise. As a result from this, slower readout speeds are used in order to not add significantly more noise.




Figure 2.5.: **Readout noise.** The pixels are shifted row-wise into the readout register, depending on the vertical shift speed ( $v_{ss}$ ) and then moved pixel-by-pixel with the horizontal shift speed into the analog to digital converter. Since noise reduction is important, minimal horizontal shift speeds will be used, while the vertical shift speed does not seem to affect the variance. To make the readout the dominant noise source, temperature was set to  $-69^\circ\text{C}$  and exposure to 10 ms. The results have been received from taking the mean and variance of multiple sets of exposures.

### 2.2.5. Quantum efficiency

When selecting cameras for scientific imaging, one of the attributes to look out for is quantum efficiency (QE). Not all incoming photons are converted into electrons, which adds an uncertainty and additional noise to measurements. This is characterized by quantum efficiency, which is defined as the fraction of photons that are converted into electrons. A QE of 100% would mean, that every incident photon is converted into an electron on the chip. This can be put into equations as

$$QE = \frac{N_{detected}}{N_{total}}, \quad (2.3)$$

where  $N_{detected}$  describes the number of detected photons and  $N_{total}$  the number of photons that were sent from the laser beam.

The measurement of the quantum efficiency has been previously carried out by Carmen Renner [5] in her diploma thesis. In order to do that, the detected photons can be rewritten according to [8]:

$$N_{total} = \frac{E_{beam}}{E_{\gamma}} = \frac{Pt_{exp}}{E_{\gamma}} \quad (2.4)$$

where the energy of all photons reaching the camera  $E_{beam}$  is described with the power  $P$  and exposure time  $t_{exp}$ . The photon energy  $E_{\gamma}$  is calculated from  $E_{\gamma} = h\nu$ .

Now the quantum efficiency reads as

$$QE = \frac{h\nu N_{detected}}{Pt_{exp}}. \quad (2.5)$$

In order to find the QE, one would now measure the accumulated counts for several exposure times, since the frequency  $\nu$  and the power  $P$  are known. For the Andor camera used in our experiment, they have been measured as 70.7 % for 671 nm and 69.9 % for 852 nm.

This camera was especially chosen from the quantum efficiency, since it is important in our case, to be able to detect most photons from Lithium and Caesium absorption.

## 2.3. Mechanical shutter

Due to the high quantum efficiency, the camera is very sensitive to stray light. Therefore it is necessary to cover the chip from stray light. A mechanical shutter which can be electronically controlled is therefore built in front of the light path, which prevents the camera from

damage between measurements. The electronics were set up and the shutter characterized in order to find perfect timings, so that light only enters the chip, when a measurement is running.

### 2.3.1. Electronic and mechanical setup

The high QE of the Andor camera translates into a high sensitivity of the CCD chip to stray light. Therefore, in order to not unnecessarily illuminate the chip between measurements, a mechanical shutter was built into the optical path, which can also be seen in Figure 2.1. By opening the shutter shortly before the imaging sequence and closing it immediately afterwards, the stray light that is implying the CCD detector can be minimized.

cite andor  
manual  
somewhere

The shutter has five fans, as seen in [, which are mechanically guided, such that they together perform a circular motion outwards](#). The circular motion is best achieved for narrow fans, therefore needing more in order to close the shutter properly.

appendix  
image

The guides are connected to each other and can be pulled outwards with a mechanical switch, that can be manually pushed or pulled. To drive the shutter electronically, a magnetic coil and a magnet are used. When the coil receives a current, the magnet, which is connected to the switch will pull the guides, moving them outwards, therefore opening the shutter.

Although the coil can pull the magnet in, it cannot be pushed away. This is compensated using a spring, which is connected to the switch. This also means, that the current driving the coil needs to be high enough, to also work against the spring.

The most optimal case would be now to have fast opening and closing times, since again we want to prevent illumination between measurements. This can be achieved by testing several springs and to make the pull from the magnetic coils as fast as possible, which means increasing the current to drive it.

A custom circuit as outlined in Figure 2.6 is used. In order to serve a high current to increase the pull from the magnetic coil on the magnet, a transistor is used which is controlled by a MOSFET driver. The complete circuit also contains a voltage regulator, so that the user does not need to know the input of the driver and can simply serve a high voltage to open the shutter ( $V \in [2, 8] \text{ V}$ ) and a low voltage ( $V = 0 \text{ V}$ ) to close it.




Figure 2.6.: **Electronic circuit to control the shutter.** The first element in the circuit is the optocoupler which is used to decouple parts from the circuit in order to not create loops, which would induce charges when the magnetic fields are on. The remaining parts can be skipped, by setting the jumpers. The inverter will flip the sign on the voltage. The MOSFET driver is used in order to serve the correct voltage at the transistor's gate and to help during discharge, such that the MOSFET will not overheat. The transistor will finally serve a high current of 2.5 A, that is throughput to the shutter. The jumpers in this figure are set as they are used in the experiment. The complete circuit can be found in Appendix B.

### 2.3.2. Dynamical properties

The shutter operates by controlling a magnetic coil, which pulls a magnet. The magnet is reverted into its original position by a spring, closing the shutter. It was discussed before, in Subsection 2.3.1, that more fans give a better approximation to a circular motion. Optimally, we would expect the shutter to open perfectly circular with linear velocity. Therefore, to optimize on this, one would look at opening and closing speeds, which should be minimal.




Figure 2.7.: **Probing the shutter for dynamics.** The shutter is probed at several positions using a laser. The offset in  $y$  direction was varied to find opening and closing times as a function of their offset.

To find the actual dynamics of the shutter, an experiment was set up. The shutter was probed at several positions, using a laser, to find the opening time. An example of one

of many measurements is shown in Figure 2.8. The laser beam has a gaussian intensity distribution, which results in a error function on the photodiode as only partial intensity is received from the diode, when the shutter is still blocking parts of the beam.

The points have been found by pointing a laser beam at a photodiode, which was blocked by the shutter. Since the laser has a finite radius, there is a transition in the signal from the minimum to the maximum, which is due to the approximately linear opening velocity of the shutter and the gaussian intensity distribution of the laser beam.




Figure 2.8.: **Shutter characterization.** The dynamics of the shutter were measured using a laser beam with a variable horizontal offset, which is fixed in this plot, and a photo diode measuring the laser intensity. In this figure, the offset is 9.3 mm from the center. An error function was fitted yielding the time until the shutter opens. The opening time was hereby defined, where the error function has its mean value, in this example being 5.2 V.

In order to find the opening and closing times of the shutter as seen in Figure 2.9a., the offset from the centre of the shutter was varied. For larger offsets, the time from the initial trigger until a signal in the photodiode was received takes longer when the shutter is opening, while it takes shorter when it is closing. In Figure 2.9b. the speed is not perfectly linear. The deviations close to the centre originate from the shape of the fans. As there are five fans, there is no mirror symmetry, so that opening to the left and right is different in the end.

These measurements will later be used to time the triggers in the measurement, so that the shutter can stay closed as long as possible. But in order to optimize it, there were two shutters at hand, and several springs to choose from, the combination of which were all tested until the timings for opening  $t_{open} = 0.12$  s and for closing  $t_{closing} = 0.14$  s were found.




Figure 2.9.: **Sample dynamics.** The figure shows opening and closing time as well as the opening velocity for a strong spring. The spring closed the shutter quickly while it prevented it from opening fast. The velocity was measured by using the beam diameter as the distance the shutter needed to transverse. For each offset, a set of 100 images were taken and the errors found as being the variance. It is noticeable, that the opening velocity on the right side is faster at first than on the left side. This is due to the structure of the shutter, as can be seen in [Appendix image of shutter]. The overall opening speed on the other hand is not affected by this and seems to be linear with the offset.

It has been found that the opening time consists of actually three timings. The initial trigger signal runs through the electronics, until it is sent to the shutter coil. The coil then has a delay, until it starts pulling the magnet. The opening time is then found by also adding the time until the fans are fully open. This is also shown in Figure 2.10.




Figure 2.10.: **Shutter opening timings.** In the figure, the four signals show the timings until the shutter is open. The blue signal is the trigger that is sent from the user to start the opening sequence. The electronics add some minor delay, which is about a few nanoseconds until they send out the green trigger. After a certain delay time, which is most likely due to the coil needing to accumulate charge first, the red signal indicates when the shutter starts to open at the centre. The opening sequence is finished, when the shutter is fully open, which is shown by the turquoise line.

## 2.4. Mask for the CCD sensor

One of the key features of the Andor iKon M is the fast kinetics readout mode. This allows for fast acquisition, which is important when imaging multiple species. Using this acquisition mode, we are able to take images  $500\text{ }\mu\text{s}$  apart from each other, so that it is possible to take consecutive images from the same cloud before it is too dilute. The acquisition mode is explained in these sections as well as the diffractions introduced by a slit, which is needed to mask parts of the chip.

### 2.4.1. Fast kinetics mode

The fast kinetics mode allows the image acquisition timings to be only dependent on the vertical shift speed, reducing the acquisition time significantly. In this mode, only a portion of the CCD is illuminated, while the dark parts of the chip will be used as a storage. This means, that as soon as an image is taken, the illuminated pixels are shifted vertically behind a mask, such that no photons can reach them any more. When the chip is full or the user has finished their acquisition, the readout process is started.

As explained earlier in Chapter 2.2.1, the readout consists of first shifting a row into the readout register and then horizontally shifting them into the ADC. This is a very time con-

suming process, since the total readout time is described by

$$t_{ij} = iv_{ss} + (i - 1)j_{max}h_{ss} + jh_{ss}, \quad (2.6)$$

where  $v_{ss}$  and  $h_{ss}$  indicate the vertical and horizontal shift speed respectively. With this equation, the readout time until a pixel with the coordinates  $i$  and  $j$  is shifted into the readout register can be calculated. The readout of all pixels beforehand is also taken into account by  $j_{max}$  which is the width of the chip in pixels (1024 for the iKon M camera). As can be seen, in the terms of the horizontal speed, the position of the last pixel in a row has a quadratic dependency. This is the dominant contribution to the readout time, which is necessary to shift the pixels into the ADC. At this point, the speed should not matter too much any more, because the experiment is already finished and the pixels are not illuminated any more.




Figure 2.11.: **Schematics of imaging in the fast kinetics mode.** **a.** Laser photons excite electrons in pixels, creating the first absorption image. **b.** The illuminated pixels are shifted down behind a cover, while the laser is shut off. **c.** The second absorption image is taken without affecting the last measurement. **d.** Both images are shifted down before starting the readout process, to not falsify the data with stray light. **e.** The readout process starts, the first row is shifted into the readout register. **f./g.** The first pixels are shifted into the ADC. It is repeated until the complete chip has been read out.

To set up the fast kinetics mode, there are several parameters that need to be set in advance.

- **Series length.** The number of images acquired before the readout phase is initiated.



- **Exposed rows.** The height of an image in the fast kinetic series.
- **Offset from bottom.** Number of rows from the bottom of the chip, which are used as a temporary storage for the illuminated pixels.

It can be easily seen that the parameters are limited by the height of the CCD in pixels:

$$\text{Exposed rows} + \text{Offset from bottom} \leq 1024 \text{ px} \quad (2.7)$$

and

$$\text{Series length} * \text{Exposed rows} \leq \text{Offset from bottom}. \quad (2.8)$$

Therefore, in our application, we chose to expose 204 px with an offset of 820 px, which gives a series length of 4. The pixels that are used as temporary storage of accumulated charge until the readout process begins need to be shielded from light. This is achieved with a slit in front of the camera that blocks part of the light beam implying on the CCD detector.

Since we will be imaging two atomic species, an absorption image for each species will be followed by two division images, as well as background images to subtract the noise.

### 2.4.2. Frequency response of a slit

The slit is an optical element in the path, which will introduce diffraction — an effect due to the wave nature of light. Figure 2.12 sketches the systematics behind a plane wave approaching a slit.

The incoming planar wave can be described by the formula:

$$E(z) = E_0 e^{-ikz}, \quad (2.9)$$

where the wave is propagating in the  $z$ -direction with the wave number  $k$  and an amplitude  $E_0$ .

It is known from Huygens principle, that each point from a planar wave can be seen as the origin of a spherical wave. The spherical waves will interfere with each other and because of the superposition principle a new wavefront will be built up.




Figure 2.12.: **Diffraction on a single slit.** A planar wave with the wavelength  $\lambda$  approaches a slit with a width  $a$ , where the maxima of the wave are drawn as dashed lines. After the slit, Huygens principle is sketched as spherical waves, (purple and green), interfering with each other. The new wavefront is created where the waves cross each other and is as such visible on the screen to the right.

The mathematical formalism of Huygens principle is simply the sum of all spherical waves, which for infinitely small distances is given as [9]:

$$E(x, z) = E_0 C \int_{Slit} \frac{1}{r^2} \exp(-ikr) dx'. \quad (2.10)$$

Here,  $r = \sqrt{(x - x')^2 + z^2}$  is the radius of a two dimensional wave,  $C$  is a normalization constant. Taylor expansion of the radius for  $z$  allows us, to substitute  $r^2$  with  $z^2$  and  $r$  with  $\frac{(x - x')^2}{2z} + z$ , while also the approximation  $(x - x') \ll z$  needs to be fulfilled. Therefore we get:

$$E(x, z) = E_0 e^{-ikz} \sqrt{\frac{ik}{2\pi z}} \int_{Slit} \exp\left(-\frac{ik}{2z}(x - x')^2\right) dx'. \quad (2.11)$$

The normalization was derived from the gaussian normal distribution  $\exp(-\frac{ik}{2z}(x - x')^2)$ . The integral will run over the slit size, with the origin in its middle and width  $a$ . We want to also write the exponential as a function of  $-\frac{i\pi t^2}{2}$ , due to the definition of fresnel integrals. The substitution follows:

$$\frac{k}{2z}(x - x')^2 = \frac{\pi t^2}{2}, \quad (2.12)$$

$$\Rightarrow t = \sqrt{\frac{k}{z\pi}}(x - x'), \quad (2.13)$$

$$\Rightarrow dx' = -dt \sqrt{\frac{z\pi}{k}}, \quad (2.14)$$

such that the field becomes

$$E(x, z) = -E_0 e^{-ikz} \sqrt{\frac{i}{2}} \int_{t(-a/2)}^{t(a/2)} \exp\left(-\frac{i\pi t^2}{2}\right) dt. \quad (2.15)$$

The fresnel integrals are defined by

$$C(x) = \int_0^x \cos\left(\frac{\pi t^2}{2}\right) dt, \quad (2.16)$$

$$S(x) = \int_0^x \sin\left(\frac{\pi t^2}{2}\right) dt, \quad (2.17)$$

so that in combination with Eulers equations and splitting up the integrals the equation for the electric field is:

$$E(x, z) = E_0 \sqrt{\frac{i}{2}} e^{-ikz} [C(-a/2) - C(a/2) - iS(-a/2) + iS(a/2)]. \quad (2.18)$$

In a real experiment, the light on a CCD chip is detected as intensity, which is given by

$$I(x, z) = \frac{2|E(x, z)|^2}{\epsilon_0 c}, \quad (2.19)$$

with the dielectric constant  $\epsilon_0$  and the speed of light  $c$ .

To verify the theory, we set up a simple experiment with a collimated laser beam pointing at the CCD of the Andor camera, with a slit in between. Such behaviour was also observed in our measurements. Figure 2.13 shows the diffraction pattern that was recorded with the Andor camera and a slit. A large, collimated beam was used to illuminate a slit in front of the CCD detector. Well resolved interference fringes can be observed behind the slit.

The experimental diffraction pattern matches the theory very well, since the deviations are in the expected regimes. The function washes out as it approaches the centre of the chip. This is due to the nature of the pixels, which only have a finite size and the fact, that the frequencies of the oscillations are lower on the outer ends. All oscillations that fit into one pixel are averaged, therefore diffraction is not visible.

This result helps to minimize the diffraction on the chip and serves as a starting point in order to optimize this effect.




Figure 2.13.: **Measuring diffraction on a slit.** In order to characterize the diffraction by the slit, the CCD detector was placed as close as possible. The parameters were measured using a ruler and yielded distance  $d = (10.9 \pm 0.5)$  mm, opening  $a = (2.5 \pm 0.5)$  mm. The wavelength was  $\lambda = 852$  nm found from the laser specifications. The blue curve is the experimental data, while the red curve was fitted, leaving the distance and opening as free fitting parameters. They were found to be  $d' = (11.0 \pm 0.3)$  mm and  $a' = (2.470 \pm 0.013)$  mm, which is in close agreement with the measured values. The residual deviation in the amplitude are caused by Gaussian intensity distribution of the laser beam used for illumination.

### 2.4.3. Optimization of the masking setup

An important issue in the optimal placement of the slit is the distance from the CCD detector at which it is going to be placed. The frequency of the diffraction pattern depends on the distance  $d$ ; it is larger for smaller  $d$ , and smaller for larger  $d$ . In the limit  $\frac{a^2}{d\lambda} \ll 1$  the diffraction pattern of a point source is recovered (Fraunhofer diffraction in the far field).

This dependence on the distance was also tested experimentally and the results are shown in Figure 2.14. The signal in one pixel is the average over all oscillations, that fit into the width of the pixel. Therefore the optimal slit position is close to the chip, since then the frequencies are maximal.




Figure 2.14.: **Distance dependent diffraction.** **a.** A slit was placed on a movable platform and diffraction was measured for various offsets  $z$ , while the slit opening  $a$  was kept constant. **b.** The diffraction frequency rises as the distance gets closer to the CCD.



### 3. Thermometry of an ultracold ideal fermi gas

The purpose of the camera is to measure scientifically important data from dense atomic clouds consisting of Lithium and Caesium. The improvement of the resolution in the whole imaging setup now allows to explore new attributes that could not be measured before and as an example, ultracold ideal fermi gases were chosen. They are apparent at very low temperatures, where only few atoms due to evaporative cooling are present. Their density distribution differs slightly from a gaussian form and from that one can extract temperatures and atom numbers.

In order to image the atoms, a technique called absorption imaging is used which is explained in the next section, before the introduction to ideal fermi gases.

#### 3.1. Absorption imaging

In order to find microscopic attributes of atoms, or systems of atoms, it is necessary to look at the atoms themselves. This is commonly accomplished using either fluorescence or absorption imaging [8]. In both cases, a laser beam is pointed at an atomic cloud, that is cooled and confined in a trap. In fluorescence imaging, the scattered light is collected, typically in a direction that is different than the illuminating beam. The intensity from the light through this method is not very high, since it is radiated in all directions. Therefore, long exposure times are required during which atoms can move and the information about the initial density and energy distribution is lost. Nevertheless, it is useful for single- and few-atom detection.

In contrast, in absorption imaging [10], the transmitted intensity of the imaging beam is recorded. Without atoms, one would see a beam profile of the laser beam. With atoms, a shadow is visible due to the atoms "blocking" the light. This is accomplished, by correctly tuning the laser to a resonance frequency of the atoms, which enables them to absorb the

cite phd thesis of Serwane, Selim, Science 332, 336 (2011)

light, exciting them to a higher state. Through spontaneous emission, the atoms will decay, making it possible to excite them once again. This method works well, when the "signal" from the absorbed light is significantly larger compared to the noise sources.

example two  
level atom  
excitation  
image

There are a set of optical elements in the imaging path, like lenses to collimate the image and refocus it, or mirrors to guide the light into the camera. Since the surfaces will most likely introduce errors into the imaging, for example from impurities or dust, only the absorption image will not suffice to gain reliable data. This is compensated by taking a total of three pictures, in order to extract only the relevant information from the image.

This can be understood when looking at the light intensity  $I_{CCD}$  reaching the camera. The atom cloud has an optical density  $OD$ , therefore the intensity can be written as

$$I_{CCD} = I_0 e^{-OD} + I_{back}, \quad (3.1)$$

where it decreases from the incident laser intensity  $I_0$  due to light scattering by atoms. The intensity  $I_{back}$  describes the background signal, that is found when the CCD is not being illuminated by a laser such as readout noise, dark noise or stray photon light. All the interesting attributes of atoms are found by looking at the optical density, therefore in order to extract that, a background frame is subtracted from the absorption image and the laser profile divided, leaving

$$\frac{I_{CCD} - I_{back}}{I_0} = e^{-OD}. \quad (3.2)$$

The laser intensity  $I_0$  is measured in a separate frame, containing the laser intensity  $I'_0 = I_0 + I_{back}$  and also the background  $I_{back}$ . Finally, the equation yields

$$\frac{I_{CCD} - I_{back}}{I'_0 - I_{back}} = e^{-OD}. \quad (3.3)$$

three images  
as example

read com-  
ment

From the resulting optical density, one can now conclude, for example, atom density distributions, atom numbers or excitation rates.



## 3.2. Density distributions of ideal Fermi gases

Ideal Fermi gases offer a new aggregate, that is complementary to Bose-Einstein condensate. They are found on the Bardeen-Cooper-Schrieffer side of the Feshbach resonances as a polarized species with only one spin component. At very cold temperatures, they start to differ from a gaussian distribution which is further investigated in this chapter.

The distribution of the atoms depends on the fraction of their temperature to the Fermi temperature [11]  $T_F$ . For  $T/T_F \gg 1$ , the atoms will follow a gaussian distribution, which can be identified using the gaussian radius:

$$\sigma_i = \sqrt{\frac{2k_B T}{m w_i^2}}, \quad (3.4)$$

with the mass of Lithium  $m$  and the trapping frequency  $w_i$  in the direction of the radius. Due to the alignment of the dipole trap, the cloud will not have a spherical shape and will therefore have different radii.

In the degenerate regime however, for  $T/T_F \ll 1$ , the radius is described by the Fermi radius

$$R_{Fi} = \sqrt{\frac{2E_F}{m w_i^2}}, \quad (3.5)$$

using the Fermi energy  $E_F$ .

It is therefore suggested to use a unified radius, as the temperatures are not known a priori:

$$R_i^2 = \frac{2k_B T}{m w_i^2} f(e^{\frac{\mu}{k_B T}}). \quad (3.6)$$

The interpolation function  $f(x)$  is hereby:

$$f(x) = \frac{Li_1(-x)}{Li_0(-x)} \quad (3.7)$$

where  $Li_n$  is polylogarithm and can be defined as

$$Li_s(z) = \sum_{k=1}^{\infty} \frac{z^k}{k^s}. \quad (3.8)$$

In our case, we integrate over all but one axes, and therefore find the fitting function for

the atom numbers:

$$n_{1D}(x) = n_{1D,0} \frac{Li_{5/2} \left( \pm \exp \left[ q - \frac{x^2}{R_x^2} f(e^q) \right] \right)}{Li_{5/2}(\pm e^q)}. \quad (3.9)$$

The derivation can be found in [11]. The parameter  $q = \frac{\mu}{k_B T}$  can be extracted, which contains information about the chemical potential  $\mu$  and the temperature  $T$ .

This parameter can then be used to calculate the degeneracy parameter:

$$\frac{T}{T_F} = [-6Li_3(-e^q)]^{-1/3}. \quad (3.10)$$

To compensate for the finite resolution of the chip on the camera, the cloud can be imaged after a short time of flight  $t$ . The temperature can then be calculated from the dynamics as

$$k_B T = \frac{1}{2} m w_i^2 \frac{R_i^2}{1 + w_i^2 t^2} \frac{1}{f(e^q)}. \quad (3.11)$$

### 3.3. Finding properties of the Fermi gas

As seen before, from an ideal Fermi gas, one can deduce the temperature and Fermi temperature of a gas. In order to implement ideal Fermi gases, Lithium was prepared in an optical dipole trap. The power was ramped down several times at the Feshbach resonance  $B = 896$  G, until only a fraction of the atoms remained. The spin-down atoms received a short laser pulse, so that the cloud only consisted of spin-up  ${}^6\text{Li}$ .

To compare the properties of the gas, the dipole trap was ramped down to three different powers, and everything I would write from here on is probably not correct, so I will be discussing this in person tomorrow as my results are super weird.




Figure 3.1.: **saywhat.** plt pls.



## **4. Conclusion and outlook**



## A. Acquisition sequence

During absorption imaging, the Andor camera will take eight images, which consist of two absorption, two division and four background images (or if the user specifies to only image one species, then there is one absorption and division image and six background images).

As the camera is only interested in the rising edge of the signal, the trigger length is not really important. But the chip will be exposed for the duration of the exposure time, after which the laser should be turned off. The, in our case, 204 pixels are then shifted downwards, taking  $t_{vshift}$ , which can be calculated by multiplying the vertical shift speed with the pixel height. For the fastest shift speed, this would result in

$$t_{vshift} = 2 \mu\text{spx}^{-1} * 204 \text{ px} = 408 \mu\text{s}, \quad (\text{A.1})$$

for the slowest speed we find

$$t_{vshift} = 64 \mu\text{spx}^{-1} * 204 \text{ px} = 13.1 \text{ ms} \quad (\text{A.2})$$




Figure A.1.: **Trigger signals.** The triggers are sent as a rectangular signal. The camera will interpret the rising edge of the signal as the starting point of the exposure. Here,  $t_{vshift}$  is the time it takes to shift the illuminated pixels downwards,  $t_{readout}$  is the time it takes the chip to read out the data and  $t_{exposure}$  the exposure time.

I'm not sure if this chapter belongs in the appendix. It is kind of a technical thing and I'm not sure where I could fit it in

Trigger number	Trigger time
1	0
2	$t_{trig}$
3	$2 * t_{trig}$
4	$3 * t_{trig}$
5	$4 * t_{trig} + t_{readout}$
6	$5 * t_{trig} + t_{readout}$
7	$6 * t_{trig} + t_{readout}$
8	$7 * t_{trig} + t_{readout}$

Table A.1.: **Trigger timing.** The trigger signals are limited by the readout speeds of the chip. This table lists the minimal timings necessary to fully read out the chip. Between signal 4 and 5, the illuminated pixels are still shifted downwards, before the chip is read out. This prevents introducing unnecessary errors, as all relevant pixels are then behind the cover of the slit.

Therefore the four subsequent signals can be taken after

$$t_{trig} = t_{vshift} + t_{exposure}. \quad (\text{A.3})$$

Four images are acquired until the readout process is started. This takes significantly longer, but as explained in Chapter 2.4.1 can be calculated via:

$$t_{ij} = i * v_{speed} + (i - 1) * j_{max} * h_{speed} + j * h_{speed}. \quad (\text{A.4})$$

For the whole chip, this gives us  $i = 4 * 204 = 816$ ,  $j = 1024$ ,  $j_{max} = 1024$ , therefore

$$t_{readout} = 816 * v_{speed} + 816 * 1024 * h_{speed}. \quad (\text{A.5})$$

The trigger signals can finally be calculated via Table A.1.



## B. Shutter circuit

The circuit given on the next page is the complete circuit as opposed to Figure 2.6. The key difference here, is that the actual board contains two input/output ports and there are also converters (visible in the top-left corner), that keep the voltage at 24 V, which is necessary in order to use the MOSFET driver.

The breadboard on the bottom are built into the circuit, so that custom additions can be made, although this is not used in our case.

---

Reference  
Gerhard



## C. Setup of the custom slit




Figure C.1.: **Drawing of the camera.** As can be seen in this drawing, the CCD chip is first hidden behind a cover, that also includes an internal shutter and then offset by an additional 5 mm.

The cover, which can be seen in Figure C.1 has a width of 12.5 mm, which adds additional space before the chip. The cover is mainly for a manual cap to cover the chip, when the camera is not used, and an internal shutter. Since we knew, that the internal shutter was not needed, we were able to remove the cover bringing the slit closer to the chip. Images of under the cover can be found in Appendix .

A/B/C? Im-  
age?

The holes for M4 screws were already there, so a custom plate was built on which the slit could be mounted on. In the technical drawing in Figure ??, the centre-most holes are reserved for the slit, which can be moved up and down to select the appropriate height needed for the imaging. The plate also gives the opportunity to move the whole set with the long holes in the outer-most edges.

It is also important to note, that the camera is very sensitive to stray light due to the high quantum efficiency. Therefore it is necessary to cover the laser path. The plate offers an-

other set of screw holes, which will hold a SM2-mount, therefore eliminating any gap that could allow photons to reach the camera externally.

The long path of SM2 tubes then only allows the smallest amount of stray light to enter the camera, which will have significantly lower intensity than the actual absorption image.

# List of Figures

2.1. Imaging path . . . . .	3
2.2. Schematic design of a CCD array and its readout . . . . .	5
2.3. Shifting charges in a CCD detector . . . . .	5
2.4. Dark noise accumulation on the chip . . . . .	8
2.5. Readout noise . . . . .	9
2.6. Electronic circuit to control the shutter . . . . .	12
2.7. Probing the shutter for dynamics . . . . .	12
2.8. Shutter characterization . . . . .	13
2.9. Sample dynamics . . . . .	14
2.10. Shutter opening timings . . . . .	15
2.11. Schematics of imaging in the fast kinetics mode . . . . .	16
2.12. Diffraction on a single slit . . . . .	18
2.13. Measuring diffraction on a slit . . . . .	20
2.14. Distance dependent diffraction . . . . .	21
3.1. saywhat . . . . .	27
A.1. Trigger signals . . . . .	31
C.1. Drawing of the camera . . . . .	35



# Bibliography

- [1] University of Heidelberg. <http://dc.zah.uni-heidelberg.de/lswscans/res/positions/q/form>. Accessed: 28.03.2016.
- [2] M. Repp, R. Pires, J. Ulmanis, R. Heck, E. D. Kuhnle, M. Weidemüller, and E. Tiemann. Observation of interspecies  ${}^6\text{Li}$ - ${}^{133}\text{Cs}$  Feshbach resonances. *Phys. Rev. A*, 87:010701, Jan 2013.
- [3] Chris McFee. An introduction to ccd operation. [http://www.mssl.ucl.ac.uk/www\\_detector/ccdgroup/opttheory/ccdoperation.html](http://www.mssl.ucl.ac.uk/www_detector/ccdgroup/opttheory/ccdoperation.html). Accessed: 23 March 2016.
- [4] Allied Vision. Guppy f-038 datasheet.
- [5] Carmen Renner. Design and characterization of a dual-wavelength high-resolution imaging system. Master's thesis, University of Heidelberg, 2014.
- [6] Timo Bastian Ottenstein. A new objective for high resolution imaging of Bose-Einstein condensates. Master's thesis, University of Heidelberg, 2006.
- [7] *Skript zum Versuch F36 – Teil I "Wellenfrontanalyse mit einem Shack-Hartmann-Sensor" des Fortgeschrittenen-Praktikums II der Universitaet Heidelberg fuer Physiker.*
- [8] Simon Murmann. Toward the exact detection of mesoscopic atom numbers. Master's thesis, University of Heidelberg, 2011.
- [9] Daniel A. Steck. Quantum and atom optics. Available online at <http://steck.us/teaching>, April 2014. revision 0.9.10.
- [10] Stephan Helmrich. Improving optical resolution by noise correlation analysis. Master's thesis, 2013.

- [11] Wolfgang Ketterle and Martin W. Zwierlein. Making, probing and understanding ultracold fermi gases. In W. Ketterle M. Inguscio and C. Salomon, editors, *Proceedings of the International School of Physics Enrico Fermi*, volume Course CLXIV, page 95. IOS Press, Amsterdam, 2008.